

Alkaline Phosphatase Isoenzyme Procedure

Cat. No. 5102

The Helena Alkaline Phosphatase Isoenzyme Procedure is intended for the qualitative determination of serum alkaline phosphatase isoenzymes by electrophoresis on cellulose acetate.

SUMMARY

Alkaline phosphatase (ALP) (EC 3.1.3.1.) is an enzyme which catalyzes the hydrolysis of phosphate esters at an alkaline pH. The greatest concentrations of ALP are found in bone, liver, intestine, and the placenta. However, practically every body tissue contains at least a small amount of ALP. Because of this wide distribution, limited information can be obtained from a total ALP assay. Fortunately each source of ALP produces one predominant isoenzyme and the tissue source of elevated ALP in serum can be determined by identifying the isoenzyme. The isoenzymes of ALP differ in their physicochemical and electrophoretic properties and, by taking advantage of these differences, the individual isoenzymes can be identified.¹ In addition to the liver, bone, intestinal and placental isoenzymes, other ALP isoenzymes have been identified in serum. These include fast liver (pre-liver), Regan, Nagao, PA, and renal isoenzymes.

A number of laboratory procedures have been used for the routine evaluation of the ALP isoenzymes. These include heat inactivation², inhibition with amino acids³⁻⁵, urea denaturation³⁻⁴, and electrophoresis on agarose⁶, paper⁷, starch gel⁸, polyacrylamide gel⁹ and cellulose acetate¹⁰⁻¹².

The Helena cellulose acetate method offers several distinct advantages over other identification methods. The technique is suitable for the electrophoretic fractionation of ALP in large numbers of sera, provides ease of handling, stability of supporting membranes and media, and reproducibility of results.

PRINCIPLE

The isoenzymes of alkaline phosphatase are separated according to their electrophoretic mobility on cellulose acetate in a tris-barbital-sodium barbital buffer. The colorimetric reaction occurs by hydrolyzing the indolyl dye.

REAGENTS

1. Alkaline Phosphatase Indolyl Blue Reagent (Cat. No. 5102)

Ingredients: When reconstituted as directed, the concentration of the reactive ingredients is as follows:

5-Bromo-3-Indolyl Phosphate p-Toluidine Salt	8.35 mM
2-Amino-2-methyl-1-propanol	1 M
Magnesium Chloride	1 mM
Stabilizers	

WARNING: FOR-VITRO DIAGNOSTIC USE. DO NOT INGEST.

Preparation for Use: Reconstitute each vial of reagent with 3 mL of Diluent. Mix the reagent to obtain complete dissolution. The reagent may be used as soon as reconstituted or within 48 hours.

Storage and Stability: The reagents should be stored at 2 to 8°C and are stable until the expiration date indicated on the label. The reconstituted reagent is stable 48 hours.

Signs of Deterioration: The dry, unreconstituted reagent should be uniformly off-white to light lavender in color.

2. Alkaline Phosphatase Indolyl Blue Diluent

Ingredients: 2-Amino-2-methyl-1-propanol and Magnesium Chloride

WARNING: FOR-VITRO DIAGNOSTIC USE. DO NOT INGEST.

Preparation for Use: The diluent is ready for use as packaged.

Storage and Stability: The diluent should be stored at 2 to 8°C and is stable until the expiration date indicated on the vial.

Signs of Deterioration: Discard the diluent if it shows signs of bacterial growth.

3. Electra® HR Buffer (Cat. No. 5805)

Ingredients: The buffer is a tris-barbital-sodium barbital buffer. **WARNING: FOR IN-VITRO DIAGNOSTIC USE. DO NOT INGEST.** The buffer contains barbital which, in sufficient quantity, can be toxic.

Preparation for Use: Dissolve one package in 750 mL deionized water. The buffer is ready for use when all material is completely dissolved.

Storage and Stability: The packaged buffer should be stored at 15 to 30°C and is stable until the expiration date indicated on the package. Diluted buffer is stable for two months at 15 to 30°C.

Signs of Deterioration: Discard packaged buffer if the material shows signs of dampness or discoloration. Discard diluted buffer if it becomes turbid.

4. Titan® III Plates (Cat. No. 3001)

Ingredients: Cellulose acetate

WARNING: FOR IN-VITRO DIAGNOSTIC USE

Preparation for Use: The plates are ready for use as packaged.

Storage: The plates should be stored at 15 to 30°C.

SPECIMEN COLLECTION AND HANDLING

Specimen: Serum is the specimen of choice. Plasma collected in heparin may be used. Anticoagulants containing oxalate, citrate or EDTA cannot be used because these substances inhibit the alkaline phosphatase activity.¹³

Patient Preparation: The patient should be fasting. Patients who have B or O blood group and are secretors may have an elevated ALP about two hours after a fatty meal.^{6, 12, 13, 18, 19}

Interfering Substances:

1. High concentrations of phosphate, oxalate, citrate and cyanide will inhibit ALP activity.^{13, 18}
2. Excess glycine may inhibit ALP activity by complexing magnesium.¹³
3. EDTA inhibits some of the isoenzymes of ALP. Do not use as an anticoagulant.¹³
4. Several drugs cause an enzymatic imbalance which may change the ALP level.^{13, 14}

Storage and Stability: It is preferable to refrigerate the blood specimen immediately after collection. Specimens should be separated from the red blood cells as soon as possible. It is strongly recommended that fresh serum samples be used. If storage is necessary, the serum should be stored frozen (-20°C) for no more than 24 hours.^{6, 15, 18}

PROCEDURE

Materials Provided: The following materials are necessary for use in the Alkaline Phosphatase Isoenzyme Test:

Item	Cat. No.
Super Z-12 Applicator	4090
Super Z-12 Sample Well Plate (2)	4096
Super CPK Aligning Base	4094
Titan Gel Chamber	4063
Dialomatic Microdispenser and Tubes	6210
1000 Staining Set	5122
Development Weight	5014

Bufferizer	5093
Titan® III Iso-Vis	3001
Electra® HR Buffer	5805
Alkaline Phosphatase Isoenzyme Control	5139
Alkaline Phosphatase Indolyl Blue Reagent	5102
Titan Blotter Pads	5037
Zip Prep	5090
Helena Marker	5000
Zip Zone® Chamber Wicks	5081
Glue Stick	5002
TITAN GEL Incubation Chamber	4062
Titan Plus Power Supply	1504


SUMMARY OF CONDITIONS

Plate	Titan® III Iso-Vis
Buffer	Electra® HR
Soaking Time for Plates	30 minutes
Sample Size	10 µL
Number of Applications	2-3
Electrophoresis Time	20 minutes
Voltage	180 V
Incubation Time	25 minutes
Incubation Temperature	37°C
Drying Time	10 minutes at approximately 56°C

STEP-BY-STEP METHOD


A. Preparation of Titan® III Plate

1. Properly code the required number of Titan® Plates by marking on the glossy, hard side with a Helena Marker. Place the mark in a corner of the plate.
2. Dissolve one bag of Electra® HR Buffer in 750 mL deionized water.
3. The plates should be soaked in the Bufferizer for 30 minutes according to the instructions for use. Alternately, the plates may be wetted by slowly and uniformly lowering a rack of plates into the buffer. The same soaking buffer may be used for soaking up to 12 plates, or for approximately one week if stored tightly closed. Improper storage may cause poor separation of the isoenzymes.





B. Preparation of Zip Zone® Chamber

1. Pour approximately 100 mL of buffer into each of the outer sections of the chamber.
2. Wet two chamber wicks in the buffer and drape one over each support bridge, being sure it makes contact with the buffer and that there are no air bubbles under the wick.
3. Cover the chamber to prevent buffer evaporation. Discard the buffer after use.





C. Sample Application

1. Place 10 µL of sample into each well of the Sample Well Plate using the Microdispenser. Cover the Sample Well Plate with a glass slide if the samples are not used within 2 minutes.
2. Prime the Super Z Applicator by quickly depressing the tips into the sample wells 3 or 4 times. Apply this loading to a piece of blotter paper. Do not load the applicator again at this point,



but proceed quickly to the next step. Priming the applicator makes the second loading more uniform.

3. Remove the wetted Titan® III Plate from the buffer with the fingertips and blot once firmly. Place a drop of water on the center of the aligning base to prevent the plate from shifting during the superimposed applications. Place the plate in the Aligning Base, cellulose acetate side up, aligning the edge of the plate with the black scribe line marked "CENTER APPLICATION". The plate should be positioned so that the identification mark is always aligned with sample #1.
4. Apply the sample to the plate by depressing the applicator tips into the sample wells 3 or 4 times and promptly transferring the applicator to the Aligning Base. Press the button down and hold it 5 seconds. Make 2 to 3 superimposed applications by repeating this step.



D. Electrophoresis of Sample Plate

1. Quickly place the plate in the chamber cellulose acetate side down. Place a weight (glass slide, etc.) on the plate to insure contact with the wicks.
2. Electrophorese for 20 minutes at 180 volts.

E. Visualization of Isoenzyme Bands

1. Remove the sample plate from the chamber at the end of the electrophoresis period and blot lightly. Place the plate, cellulose acetate side up, on the blotter. Pipette 1.5 mL of the reagent onto the cellulose acetate surface. Tilt the blotter until the reagent covers the surface of the plate. Allow the reagent to soak into the plate for 1 minute.
2. Lay a clean glass rod or a serological pipette on the cellulose acetate, and gently roll it across the plate to remove any excess reagent. Failure to remove sufficient excess reagent or excessive pressure on the rod will cause smearing of the pattern.
3. Place the plate, acetate side up, into a preheated Incubation Chamber for 30 minutes at 37°C.
4. After incubation, place the plate in a staining rack, and immerse it in 5% acetic acid for 5 minutes.
5. Then, immerse the plate and rack in water for 5 minutes.
6. Remove the plate from the rack and lay it on a blotter. Dry it in a 56°C oven for 10 minutes.

F. Evaluation of the ALP bands

Qualitative evaluation: The ALP plates may be inspected visually for the presence of the isoenzyme bands.

Stability of end product: The plates should be viewed within the working day. Protect from light in the interim.

Quality Control: The Alkaline Phosphatase Isoenzyme Control (Cat. No. 5139) verifies all phases of the procedure and should be used on each plate run. The control may be used as a marker for the proper location of the bands.

RESULTS

Results of the alkaline phosphatase isoenzyme migration are described in comparison to typical serum protein migrations. Fast liver (pre-liver) migrates in the protein alpha₁ region and the major liver band migrates in the alpha₂ region. Placental, Regan, Nagao, and renal isoenzymes migrate with bone in the alpha₂/beta region, but they appear as tighter bands than bone. PA migrates cathodic to the intestinal band in the gamma region. An ultra-fast band which migrates in the albumin region is seen occasionally using the visualization method. It is believed to be caused by bilirubin bound to albumin.^{29, 30} A specimen containing both liver and bone may exhibit one wide diffuse band in the alpha₂ pre-beta region. When such a pattern is obtained, heat inactivation may be helpful in differentiating the two isoenzymes.

true isoenzyme. Further research must be conducted to determine its true origin and significance.

BIBLIOGRAPHY

1. Fishman, W.H., Am J Med, 56:617-650, 1974.
2. Posen, S. et al., Ann Int Med, 62(6):1234-1243, 1965.
3. Stepan, J. and Vecerek, B., Acta Univ Carol, Mono 77:135-140, 1977.
4. O'Carroll, D. et al., AJCP, 63:564-572, 1975.
5. Rhone, D.P. and Mizuno, F.M., AJCP, 59:531-541, 1973.
6. Sundblad, L. et al., Clin Chim Acta, 45:219-223, 1973.
7. Baker, R.W.R. and Pellegrino, C., Scan J Clin Lab Invest, 6:94-101, 1954.
8. Keiding, N.R., Scan J Clin Lab Invest, 11:106-112, 1959.
9. Johnson, R.B., et al., Clin Chem, 18(2):110-115, 1972.
10. Fritsche, H.A. and Adams-Park, H.R., Clin Chem, 18(5):417-421, 1972.
11. Viot, M. et al., Biomedicine, 31:74-77, 1979.
12. Rhone, D.P. et al., Clin Chem, 19(10):1142-1147, 1973.
13. Young, D.S. et al., Clin Chem, 21(5):246D-248D, 1975.
14. Ahmad, I., Clin Chem, 24(10):1850-1851, 1978.
15. Massion, C.G. and Frankenfeld, J.K., Clin Chem, 18(4):366-373, 1972.
16. Dysert, P., Personal Communication, May, 1981.
17. Lee, L.M. and Kenny, M.A., Clin Chem, 21(8):1128-1135, 1975.
18. Wolf, P.L., Arch Pathol Lab Med, 102:497-501, 1978.
19. Eastman, J.R., and Bixler, D., Clin Chem, 23(9):1769-1770, 1977. 20. Cherian, A.G., and Hill, J.G., AJCP, 70(5):783-789, 1978.
21. Nerenberg, S.T., Medical Technology, Lea and Febiger, Inc. Philadelphia 122-123, 1973.
22. Guilleux, F., et al., Clin Chim Acta, 87:383-386, 1978.
23. Fishman, W.H., et al., AJCP, 57:65-74, 1972.
24. Messer, R.H., Am J Obst & Gynec, 98(4):459-465, 1967.
25. Kranc, M.C. et al., Clin Chem, 24(6):1012, 1978.
26. Fishman, W.H. et al., Nature, 219:697-699, 1968.
27. Inglis, N.R. et al., Clin Chim Acta, 33:287-292, 1971.
28. Cha, C.M., et al., Clin Chem, 21(8):1067-1071, 1975.
29. Tsung, S.H., Clin Chem, 24(11):2068, 1978.
30. Hardin, E., et al., Clin Chem, 24(1):178-179, 1978.

For Sales, Technical and Order Information and Service Assistance, call 800-231-5663 toll free.

Helena Laboratories warrants its products to meet our published specifications and to be free from defects in materials and workmanship. Helena's liability under this contract or otherwise shall be limited to replacement or refund of any amount not to exceed the purchase price attributable to the goods as to which such claim is made. These alternatives shall be buyer's exclusive remedies.

In no case will Helena Laboratories be liable for consequential damages even if Helena has been advised as to the possibility of such damages.

The foregoing warranties are in lieu of all warranties expressed or implied including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

